




Can AI Models See What Anatomists See?

*Identifying Deep Learning Architectures Optimized for **Oral Ultrasonography**:
A Comparative Study in **Buccal Mucosa Classification***

Goun Jung, Yu-Ran Heo, Hee-Jin Kim

Division in Anatomy and Developmental Biology, Department of Oral Biology
Human Identification Research Institute, BK21 FOUR Project
Yonsei University College of Dentistry


Opening 4 Keywords




*Infraoral
Anatomy*



Ultrasonography



*Artificial
Intelligence*




Explainability

"Can deep learning models classify upper and lower buccal mucosa regions from ultrasound images in a way that aligns with anatomical understanding?"

Background

Anatomical Significance of the Buccal Mucosa

- Forms the **Buccinator Musculo-Mucosal Flap (BMMF, Bozola Flap)** with the buccinator muscle
- First choice for the reconstruction of defects in oral cavity, oropharynx, and nasal septum
- The **BMMF** is an extremely versatile 'like for like' local flap option due to its long arc of rotation
- Receives blood supply from fascial artery and buccal artery: Classified as an axial flap: classified as an axial flap




a long arc of rotation of BMMF: 2013, K.Khan et al.

Background

Clinical Advantages of the Buccal Mucosa

- Provides mucosal coverage, not cutaneous, and also maintains sensation
- A valuable reconstructive option that can cover extensive areas while maintaining vascular supply
- Donor site can be closed primarily with minimal deformity or scarring





Nine-Month Post-surgery: Palatal Mucoepidermoid Carcinoma Site
Reconstruction with BMMF - Surgical and Donor Site Healing:
2021, Sesenna et al.

Background

Preoperative Challenge


- Precise identification of **anatomical structures like arteries or duct** is critical
- **Ultrasound is a non-invasive tool**, but interpretation is difficult
 - Anatomically ambiguous region
 - Few previous studies
 - High subjectivity in image reading

An objective, AI-assisted approach is needed to assess vascular anatomy and flap viability.




presents unique challenges due to complex tissue structures and low contrast in the oral region: 2018, Yiqun Liu et al.

Research Objective



Automated Classification

To classify buccal mucosa ultrasound images into superior and inferior regions using deep learning



Visual Explainability

To visualize model attention using Grad-CAM and evaluate anatomical focus



Model Optimization

To improve model robustness and generalization via hyperparameter tuning and patient-level data separation



Anatomical Validation

To validate model predictions against known anatomical structures, ensuring clinical interpretability

Methods Workflow Overview








Image Set Description

477 intraoral ultrasound images
collected from 40 volunteers






✓ Less than 1% difference between classes






- **Subject-independent evaluation**
No patient overlap across training, validation, and test sets.
To ensure generalizability and prevent data leakage.

*Obtained IRB Approval (2024-2-0006)

Classification Results




Classification Results







Grad-CAM Results

Representative Attention Maps




IO3_IBM3_R







AvgCAM_inferior

Grad-CAM Results

Representative Attention Maps










V6_SBM3_L



AvgCAM_superior

Grad-CAM Results


Misguided Model Attention



Conclusion



Deep learning model successfully classified buccal mucosa ultrasound images



Grad-CAM intensity patterns offer insight into model behavior and error characteristics.



Potential applications in various clinical domains with Detection of key anatomical structures

" By aligning AI attention with anatomical expertise, we demonstrate how explainable deep learning can advance ultrasound interpretation and clinical decision-making."

You might ask...

? *Why CNN Models? Why ResNet-50 over DenseNet-121?*

? *No augmentation. Why?*


? *Why Grad-CAM?*

? *Ultrasound-Specific Layer?*

? *Overfitting Problem?*

? *What is Optuna?*

? *Why F1?*



 gounjung.me

- *In the age of AI, curiosity is our most powerful tool.*
- *Data is fuel, AI is the engine— but direction comes from us.*



You might ask...

? *Why CNN Models? Why ResNet-50 over DenseNet-121?*

? *No augmentation. Why?*

? *Why Grad-CAM?*

? *Ultrasound-Specific Layer?*

? *Overfitting Problem?*

? *What is Optuna?*

? *Why F1?*

QnA

? Why CNN Models? Why ResNet-50 over DenseNet-121?

We chose CNN models as they have proven effective in medical image classification, especially where anatomical structures need to be spatially interpreted.

Among the tested architectures, **ResNet-50** provided a good balance between **model complexity, training time, and interpretability**, which was critical for our Grad-CAM-based explainability.

Although **DenseNet-121** showed slightly higher accuracy, ResNet-50 was chosen for final Grad-CAM visualization because:

- Its skip connections make attention maps more localized and interpretable.
- It had more **stable convergence** during cross-validation.
- And it is **widely adopted in clinical imaging research**, supporting reproducibility.

? No augmentation. Why?


We deliberately **did not apply data augmentation** in this experiment because:

- The dataset is **anatomically sensitive**. Flipping, rotation, or shifting could disrupt spatial orientation (e.g., left vs. right cheek, superior vs. inferior).
- Our **primary goal was model interpretability**, not just accuracy.
- We focused on **inherent spatial patterns**, as learned from unaltered, real patient data.

In future work, domain-specific augmentation strategies could be explored—such as **intensity variation or speckle simulation**.


 gounjung.me

- *In the age of AI, curiosity is our most powerful tool.*
- *Data is fuel, AI is the engine— but direction comes from us.*



Why these 4 Models?

<https://medium.com/@enrico.randellini/image-classification-resnet-vs-efficientnet-vs-efficientnet-v2-vs-compact-convolutional-c205838bbf49>



◆ ResNet-50

is a well-established baseline model: moderately sized and widely used in medical imaging.

◆ DenseNet-121

offers improved accuracy with better feature reuse, though at higher computational cost.

◆ EfficientNet-B0

was selected for its **small size and lightweight architecture** : useful in clinical deployment scenarios.

◆ EfficientNet-B7

on the other hand, represents a **high-end accuracy-focused model**, allowing us to test the upper bound of performance.

Hyperparameter Optimization Results

Deep Learning Hyperparameters

Deep Learning Model

